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Measurement of number concentration of ultrafine bubbles using laser diffraction method
and comparison with particle tracking analysis method

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Abstract For promoting industrial applications of ultrafine bubbles, it is essential to establish the evaluation method of the characteristics quantitatively. We have measured the size and number concentration of ultrafine bubbles generated by GALF system using five types of ten measuring instruments from around 2009. Although number concentration measurement using laser diffraction method has been challenging, recent improvements in technique have extended its usefulness. In this study we measured the number concentration of ultrafine bubbles using laser diffraction method and compared the result with a measurement using particle tracking analysis method, performed on the same samples. We found that the number concentration of ultrafine bubbles was quite similar as determined by either method. Therefore, it was confirmed that laser diffraction method can be used for measuring number concentration of ultrafine bubbles.

Keywords: Ultrafine bubbles, Laser diffraction method, Particle tracking analysis method, ultrafineGALF, Number concentration

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Table 1 Comparison of measurement methods of fine bubbles

	Measurement methods	
	Laser diffraction	Particle tracking analysis
Measurable size (nm)	80~100,000	30~1000
Measurable number concentration (/mL)	$10^5 \sim 10^{10}$	$10^7 \sim 10^9$
Quantity of sample (mL)	7	0.5
Necessary pretreatment	No need	No need

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
Fig. 1 Schematic image of laser diffraction measurement system for fine bubbles.


Fig. 2 Data processing of quantitative laser diffraction method about fine bubbles measurement.

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
Fig. 3 Number concentration distribution of UFBs measured by particle tracking analysis method.


Fig. 4 Number concentration distribution of UFBs measured by particle tracking analysis method.

Table 2 Total number concentration of UFBs measured by laser diffraction method and particle tracking analysis method

	Sample A (E8/mL)	Sample B (E8/mL)	Sample C (E8/mL)
Laser diffraction (SALD-7500X10)	26.7	9.7	5.5
Particle tracking analysis (NanoSight NS500)	23.5	16.7	12.2

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